

(6) ~~ob~~ wait on \$

-B res p a b i l i d a d l u s a r i o e s t a b i m e r e l l a p a l i b r a t i o n e l e q p B t i e m p l i d a z a l o s r e s u l t a d o s i n f o r m a c e n e s t e d o m i n e b a s e r a t e r i s t i a s p o e q p d l a s m i t o r e s d e r a t o r y d l a s b u s p a t i a s d u s o y i a d  
( I t i s i t e p o n s b i l i t y o f t h e u s e o f t h a t i n i e b a t i o n d a t a o f h i g h e q u i p m e n t n t d e i t a d a l i y o f t e m l s t p d r i n i t i s p a c h s o n i t e k a d s o f t e q u i p m e n t o p a t n g a n d i t o n : a n d g o o d u s a g e p a i n t

**Eusoides resulta de la alibación del usuario**  
(The user of alibation is the result)

Is resultats y niveles d'inertidm bresenbatare d'itifa de alibracón res pñeizaxa ment e al instrumet oescit oenla hja 1.  
(The island itele lo un nite dilem itis adibaton ofitaton d'u nvey d'tei nsmatidhe d'attemonit adibaton)

~~-ls resulta de qe se present a nenest e eereinfoda zabilid de p t r ones nacionales.~~  
(The ~~en~~ l st p o d i n l t i s i n f a h a v e a t a b i l i t y d n a t o n a l s a n d a d s

Inertia brepne se expresa unátemente a  $k = 2$ , que se representa en la figura 9.5m%  
(The grand unit is not with accuracy of  $k = 2$ , which is represented in figure 9.5m%)

**Inertion breed medd fe estimad Neph H140 MN 2002 Cia prala exesi dition breenas medcos:**  
(Theurent oftemagnet wisiat auding dte N X - CH40IM 2002 Cia dition ofurent in Itemagnet)

(Despite on of method)

**- La biblioteca consiste en la medicina leída para la muerte del hombre es la del patrón**  
(The alibaton on istoffte editenag en bftedi sl an toffte B. vilt ite ad onff)

-la reabilitabil de a determina o intervală de medii de corelație mai mică decât cea a celor 3 indicații selectate.


[illegible]

~~The equipment shall consist of one system at all times representing 1 inch (symbol) = 1 inch (at station)~~

- ~~Geometrische~~ ~~libra~~ ~~ensistema~~ ~~inglês~~ ~~ment~~ ~~a~~ ~~els~~ ~~gient~~ ~~at~~ ~~de~~ ~~over~~ ~~sa~~ ~~in~~ ~~25,4~~ ~~mm~~ .  
(When the equipment is installed in accordance with the following conditions, the offset = 25.4 mm)

Mess Servicios Metrológicos S. de R.L. de C.V. Acceso III, No. 16A, Nave 10, Parque Industrial Benito Juárez, Querétaro, Qro.  
Para cualquier duda, comentario, sugerencia, felicitación o queja favor de llamar o comunicarse a los siguientes correos:  
Tel. (442) 1 96 49 38,    oscar@mess.com.mx    marypaz.cruz@mess.com.mx    calidad@mess.com.mx

CERTIFICADO DE CALIBRACIÓN  
DIGITAL ORIGINAL  
MESS SERVICIOS METROLÓGICOS



(Laboratory accredited by ema with accreditation number D-97 as of 2010-08-18.  
In compliance with ISO/IEC 17025:(valid) NMX-EC-17025-IMNC-valid.  
"General requirements for the competence of testing and calibration laboratories").

**"La emisión de certificados de calibración/informes de medición o ensayo es de manera electrónica en cumplimiento de los estándares"**

